

大腸菌由来の再構成型無細胞系内での翻訳反応は 大腸菌細胞と同一か？


金森 崇、布施(村上) 朋重、松本 令奈 (ジーンフロンティア(株))

2024. 6. 13

PURE system

PURE $frex$ is based on the PURE system technology.

The PURE system is a reconstituted cell-free protein synthesis system, which consists only of purified factors necessary for transcription, translation and energy regeneration in *Escherichia coli*.



Shimizu Y. *et al.* (2001) *Nat. Biotechnol.*, vol. 19, p. 751-755.
 Shimizu Y. *et al.* (2005) *Methods*, vol. 36, p. 299-304.

Topics

- ★ PURE system (Translation in *E. coli*)
- ★ Initiation in the PURE system
- ★ Optimum sequence for the PURE system
 - ▶ 5'-UTR
 - ▶ N-terminal region in the ORF


Topics

- ★ **PURE system (Translation in *E. coli*)**
- ★ Initiation in the PURE system
- ★ Optimum sequence for the PURE system
 - ▶ 5'-UTR
 - ▶ N-terminal region in the ORF

PURE system

PURE $frex$ is based on the PURE system technology.


The PURE system is a reconstituted cell-free protein synthesis system, which consists only of purified factors necessary for transcription, translation and energy regeneration in *Escherichia coli*.



Shimizu Y. *et al.* (2001) *Nat. Biotechnol.*, vol. 19, p. 751-755.

Shimizu Y. *et al.* (2005) *Methods*, vol. 36, p. 299-304.

Initiation



MTF (Methionyl-tRNA transformylase)

- transfers formyl group from 10-formyl-tetrahydrofolic acid (10-CHO-THF) to the free amino group in methionine bound to initiator tRNA

IF1 (Initiation Factor 1)

- binds the A-site in the 30S subunit and prevents binding of aminoacyl-tRNA.
- modulates and stabilizes binding of IF2 and IF3 to the 30S subunit.

IF2 (Initiation Factor 2)


- is a GTPase and hydrolyzes GTP to GDP during formation of the 70S initiation complex.
- binds formylmethionyl-initiator tRNA and promotes its entry to the P-site.
- is expressed as isoforms (α , β , β') from in-frame different AUG codons of *infB* gene.

IF3 (Initiation Factor 3)

- binds to the 30S subunit and prevents association of the premature 50S subunit.
- enhances the accommodation of formylmethionyl-initiator tRNA to the P-site.
- Its gene contains the non-canonical initiation codon AUU.

All of four proteins are essential for cell growth.

Elongation



EF-Tu (Elongation Factor Tu)

- delivers aminoacyl-tRNA to the A-site of ribosome in a GTP-bound form.
- is encoded by two genes, *tufA* and *tufB*.

EF-Ts (Elongation Factor Ts)


- associates with GDP-bound EF-Tu and stimulates the exchange of GDP and GTP.

EF-G (Elongation Factor G)

- facilitates the translocation of the ribosome in a GTP hydrolysis-dependent manner.

All of three EFs are essential for cell growth.

Elongation



EF-Tu (Elongation Factor Tu)

- delivers aminoacyl-tRNA to the A-site of ribosome in a GTP-bound form.
- is encoded by two genes, *tufA* and *tufB*.

EF-Ts (Elongation Factor Ts)

- associates GDP-bound EF-Tu and stimulate the exchange of GDP and GTP.

EF-G (Elongation Factor G)

- facilitates the translocation of the ribosome in a GTP hydrolysis-dependent manner.

EF-P and ABC-F proteins

- promotes the formation of the peptide bond for the specific sequences.

EF-P: poly-Pro

YheS: SecM, **Uup:** poly-Pro, **YbiT:** K/E repeat and D/E repeat, **EttA:** D/E repeat

Ude *et al.* (2013) *Science*, 339, 82-85.

Doerfel *et al.* (2013) *Science*, 339, 85-88.

Murina *et al.* (2019) *J. Mol. Biol.*, 431, 3568-3590.

Chadani *et al.* (2024) *Nucleic Acids Res.*, gkae309.

Elongation (Effect of EF-P)

AmiB (*E. coli*)




129

VPPPPPPPPV

IL6 (mature) (*H. sapiens*)


1

MVPPGE



EF-P enhanced the synthesis of proteins with Pro-Pro sequence in the PURE system.

Termination



RF1 (Release Factor 1)

- recognizes the stop codon UAG and UAA, and releases the nascent polypeptide chain.
- is post-translationally methylated at the position Gln 235 in the conserved GGQ motif.

RF2 (Release Factor 2)


- recognizes the stop codon UGA and UAA, and releases the nascent polypeptide chain.
- is post-translationally methylated at the position Gln 252 in the conserved GGQ motif.
- Its gene contains UGA codon at the position 26 and the full-length protein is synthesized by +1 frameshift.
- Substitution of Thr at the position 246 to Ala increases the termination efficiency.

RF3 (Release Factor 3)

- is a ribosome-dependent GTPase.
- promotes the release of RF1 and RF2 from the ribosome after the release of the peptide chain.

RF1 and RF2 are essential for cell growth.

Termination



RF1 (Release Factor 1)

- recognizes the stop codon UAG and UAA, and releases the nascent polypeptide chain.
- is post-translationally methylated at the position Gln 235 in the conserved GGQ motif.

RF2 (Release Factor 2)


- recognizes the stop codon UGA and UAA, and releases the nascent polypeptide chain.
- is post-translationally methylated at the position Gln 252 in the conserved GGQ motif.
- Its gene contains UGA codon at the position 26 and the full-length protein is synthesized by +1 frameshift.
- Substitution of Thr at the position 246 to Ala increases the termination efficiency.

RF3 (Release Factor 3)

- is a ribosome-dependent GTPase.
- promotes the release of RF1 and RF2 from the ribosome after the release of the peptide chain.

RF1 and RF2 are essential for cell growth.

Recycling



RRF (Ribosome Recycling Factor)

- Is involved in ribosome recycling with EF-G after the termination process.

RRF is essential for cell growth.


E. coli-based cell-free protein synthesis systems

	Extract system	Reconstituted system		
	S30 system	PURE system (original)	PUREfrefx [®] 1.0	PUREfrefx [®] 2.0
Typical Yield (µg/mL)	100-1,000	10-200	10-200	20-1,000
Contamination				
RNase	very High	Low	very Low	very Low
LPS	very High	High	very Low	very Low
Template DNA				
Plasmid DNA	OK	OK	OK	OK
PCR product	NG	OK	OK	OK
Customization of Reagent	Difficult	Easy	Easy	Easy
Purification of His-tagged product	OK	NG	OK	OK

Shimizu Y. *et al.* (2001) *Nat. Biotechnol.*, vol. 19, p. 751-755.

Shimizu Y. *et al.* (2005) *Methods*, vol. 36, p. 299-304.


Example of protein synthesis




Topics

- ★ PURE system (Translation in *E. coli*)
- ★ **Initiation in the PURE system**
- ★ Optimum sequence for the PURE system
 - ▶ 5'-UTR
 - ▶ N-terminal region in the ORF


Initiation in the PURE system (MTF and 10-CHO-THF)



10-CHO-THF




MTF




Formylation of initiator methionine was not essential for the translation in the PURE system.

Initiation in the PURE system (IFs)




IF1, IF2, IF3





Initiation factors were not always required for the translation reaction.

Initiation in the PURE system (initiation codon)



sfGFP
 M S K G E E ..
 atgtctaaagggtgaagaa..

mCherry
 M A I I K E ..
 atggctattattaaagaa..



Initiation codons in *E. coli* (Ref; Blattner et al. (1997) Science, 277, p. 1453)

- AUG; 83%
- GUG; 14% (e.g. EF-Tu, ArgRS, HisRS, RRF)
- UUG; 3%
- AUU; 2 genes (*infC* (IF3), *pcnB*)

GTG, TTG, and ATT functioned as the initiation codon, but the synthesis efficiency was dependent on the target proteins.

Summary

Initiation in the PURE system

- ▶ All IFs and formylation of initiator methionine are not essential for the translation in the PURE system, but IF2 and IF3 enhance the translation efficiency.
- ▶ Codons used as initiation codon in *E. coli* are available for the translation in the PURE system, albeit with low efficiency.

Topics

- ★ PURE system (Translation in *E. coli*)
- ★ Initiation in the PURE system
- ★ **Optimum sequence for the PURE system**
 - ▶ **5'-UTR**
 - ▶ N-terminal region in the ORF

5' UTR of the template DNA for PURE $frex$

Currently used sequence (derived from T7 *gene 10* UTR)

GAAATTAATACGACTCACTATAGGGGAGACCACAACGGTTTCCCTCTAGAAATAATTTTGTTTAACTTTAAGAAGGAGATATAACCAATGNNNNNNNNNNNNNNNNNN



Stemloop region is necessary for efficient transcription reaction and stabilization of the transcript.

AT-rich region binds S1 protein of 30S ribosome subunit and increases translation efficiency.

SD (Shine-Dalgarno) sequence binds 3'-terminus of 16S rRNA and localizes mRNA to the start position of translation.

SD and Spacer



Long SD sequence

< 10 nt spacer

SD	Spacer
TAAGGAGGTG	TATAATATACCA
	TAATATACCA
	ATATACCA
	TACCA
	CCA

Currently used SD sequence


5-12 nt spacer

SD	Spacer
AAGGAG	TATAATATACCA
	TAATATACCA
	ATATACCA
	TACCA
	CCA

Short SD sequence


8-10 nt spacer

SD	Spacer
AAG	TATAATATACCA
	TAATATACCA
	ATATACCA
	TACCA
	CCA



currently used sequence (control)



Stemloop



Stemloop region

3 nt


6 nt



Stemloop and AT-rich



Stemloop	AT-rich	
GGGAGACCACAACGGTTTCCC	TCTAGAAATAATTTTGTTTAACTTTAAG	
GGGAGA	TCTAGAAATAATTTTGTTTAACTTTAAG	
	TCTAGAAATAATTTTGT	G
	AATAATTTTGTTTAACTTTAAG	
	AATTTTGTTTAACTTTAAG	
	TTTGTTTAACTTTAAG	
	GTTTAACTTTAAG	
	TTAACTTTAAG	
	G	
GGGAGA	TTTGTTTAACTTTAAG	
	TTT A TTTAA A TTTAAG	
	TTT A TTTAA T TTTAAG	



currently used sequence (control)

6nt-SL

AT-rich region	No	< 12 nt	> 15 nt
	×	△	○

Effect of deletion of 5' UTR of the downstream gene


(+) AT-rich/SD



(-) AT-rich




(-) SD



(-) AT-rich/SD




(-) linker




SD sequence was required, but AT-rich region was not, for translation of the downstream gene.

Effect of deletion of 5' UTR of the upstream gene

(+) AT-rich/SD (+) AT-rich




(-) AT-rich (+) AT-rich



Short SD (+) AT-rich




(-) SD (+) AT-rich




(+) AT-rich/SD (-) AT-rich




(-) AT-rich (-) AT-rich



Short SD (-) AT-rich



(-) SD (-) AT-rich



Deletion of AT-rich region or SD sequence in the upstream gene increased translation of the downstream gene only with AT-rich region in the 5' UTR.

Summary

5' UTR for the translation in the PURE system



Currently used sequence (control)

Efficiency

GAAATTAATACGACTCACTATAGGGGAGACCACAACGGTTTCCCTCTAGAAATAATTTTGTTTAACTTTAAGAAGGAGATATACCA

100%

GAAATTAATACGACTCACTATAGGG-----TCTAGAAATAATTTTGTTTAACTTTAAGAAGGAGATATACCA

50%

GAAATTAATACGACTCACTATAGGGGAGA-----TCTAGAAATAATTTTGTTTAACTTTAAGAAGGAGATATACCA

85%

GAAATTAATACGACTCACTATAGGGGAGACCACAACGGTTTCC-----GAAGGAGATATACCA

10%

GAAATTAATACGACTCACTATAGGGGAGACCACAACGGTTTCCCTCTAGAAATAATTTTGTTTAACTTTAAG-----ATATACCA



5%

6 nt for efficient transcription

essential for translation
> 12 nt for efficient translation

essential for translation
> 3 nt for efficient translation

Summary




Topics

- ★ PURE system (Translation in *E. coli*)
- ★ Initiation in the PURE system
- ★ **Optimum sequence for the PURE system**
 - ▶ 5'-UTR
 - ▶ **N-terminal region in the ORF**

Effect of the sequence on the protein synthesis

Example of protein synthesis



Synthesis efficiency is highly dependent on the target proteins.

Effect of N-terminal sequence on the protein synthesis

ARTICLE

<https://doi.org/10.1038/s41467-019-13810-1>

OPEN

A short translational ramp determines the efficiency of protein synthesis

Manasvi Verma^{1,8}, Junhong Choi^{2,3,6,8}, Kyle A. Cottrell^{1,8}, Zeno Lavagnino^{1,7}, Erica N. Thomas⁴, Slavica Pavlovic-Djuranovic¹, Pawel Szczesny⁵, David W. Piston¹, Hani S. Zaher⁴, Joseph D. Puglisi² & Sergej Djuranovic^{1*}

Translation initiation is a major rate-limiting step for protein synthesis. However, recent studies strongly suggest that the efficiency of protein synthesis is additionally regulated by multiple factors that impact the elongation phase. To assess the influence of early elongation on protein synthesis, we employed a library of more than 250,000 reporters combined with in vitro and in vivo protein expression assays. Here we report that the identity of the amino acids encoded by codons 3 to 5 impact protein yield. This effect is independent of tRNA abundance, translation initiation efficiency, or overall mRNA structure. Single-molecule measurements of translation kinetics revealed pausing of the ribosome and aborted protein synthesis on codons 4 and 5 of distinct amino acid and nucleotide compositions. Finally, introduction of preferred sequence motifs only at specific codon positions improves protein synthesis efficiency for recombinant proteins. Collectively, our data underscore the critical role of early elongation events in translational control of gene expression.

Verma *et al.* (2019) *Nat. Commun.*

Peptidyl transferase center decompaction and structural constraints during early protein elongation on the ribosome

Bin Jia¹, Tianlong Wang^{1✉} & Jean Lehmann^{2✉}



Peptide bond formation on the ribosome requires that aminoacyl-tRNAs and peptidyl-tRNAs are properly positioned on the A site and the P site of the peptidyl transferase center (PTC) so that nucleophilic attack can occur. Here we analyse some constraints associated with the induced-fit mechanism of the PTC, that promotes this positioning through a compaction around the aminoacyl ester orchestrated by U2506. The physical basis of PTC decompaction, that allows the elongated peptidyl-tRNA to free itself from that state and move to the P site of the PTC, is still unclear. From thermodynamics considerations and an analysis of published ribosome structures, the present work highlights the rationale of this mechanism, in which the free-energy released by the new peptide bond is used to kick U2506 away from the reaction center. Furthermore, we show the evidence that decompaction is impaired when the nascent peptide is not yet anchored inside the exit tunnel, which may contribute to explain why the first rounds of elongation are inefficient, an issue that has attracted much interest for about two decades. Results in this field are examined in the light of the present analysis and a physico-chemical correlation in the genetic code, which suggest that elementary constraints associated with the size of the side-chain of the amino acids penalize early elongation events.

Jia *et al.* (2021) *Sci. Rep.*

**N-terminal 3rd to 5th amino acid residues are important for efficient translation.
Larger amino acids, which are encoded by AT-rich codon, facilitated early elongation.**


N-terminal codon suitable for PURE_{frex}

N-terminal codon of Trastuzumab HC




Yield: 0.6 mg/mL

	Met	Glu	Val	Gln	Leu	Val	GC (%)	ΔG (kcal/mol)
AT	ATG	GAA	GTA	CAA	TTA	GTT	28	0.0
Major	ATG	GAA	GTG	CAG	CTG	GTG	56	-4.7
GC	ATG	GAG	GTG	CAG	CTG	GTC	61	-3.9




AT-rich codons > Major codons

GC content




mRNA secondary structure




Effect of temperature on the protein synthesis

PUREfres 2.1 (4 mM GSH)
 + template DNA
 incubation
 for 24 hours
 at 23, 30 or 37°C
 SDS-PAGE



Fuse-Murakami *et al.* (2024) *Int. J. Mol. Sci.*, 25, 5264.



Synthesis efficiency at lower temperature depended on the target proteins.

Effect of 4th amino acid on sfGFP synthesis

sfGFP

M S K **G** E E L F T G
 ATGTCTAA**GGT**GAAGAATTATTTACTGGT..

	T		C		A		G		
T	TTT	Phe	TCT	Ser	TAT	Tyr	TGT	Cys	
	TTC		TCC		TAC		TGC		
	TTA	Leu	TCA		TAA	Stop	TGA	Stop	
	TTG		TCG		TAG	Stop	TGG	Trp	
C	CTT		Leu	CCT	Pro	CAT	His	CGT	Arg
	CTC			CCC		CAC		CGC	
	CTA	CCA		CAA		Gln	CGA		
	CTG	CCG		CAG		CGG			
A	ATT	Ile	ACT	Thr	AAT	Asn	AGT	Ser	
	ATC		ACC		AAC		AGC		
	ATA		ACA		AAA	Lys	AGA	Arg	
	ATG	Met	ACG		AAG	AGG			
G	GTT	Val	GCT	Ala	GAT	Asp	GGT	Gly	
	GTC		GCC		GAC		GGC		
	GTA		GCA		GAA	Glu	GGA		
	GTG		GCG		GAG	GAG	GGG		


Effect of 4th amino acid on sfGFP synthesis

PUREfrex 2.1 (4 mM GSH)
+ template DNA
incubation
for 24 hours
at 23, 30 or 37°C

Fuse-Murakami et al. (2024) Int. J. Mol. Sci., 25, 5264.

sfGFP

M S K **G** E E L F T G
ATGTCTAAAG**GT**GAAGAATTATTTACTGGT..



The fourth amino acid was very important for the efficient translation at lower temperature.

Effect of 4th amino acid on sfGFP synthesis

PUREflex 2.1 (4 mM GSH)

+ template DNA

incubation


for 24 hours

at 23, 30 or 37°C

sfGFP

Fuse-Murakami *et al.* (2024) *Int. J. Mol. Sci.*, 25, 5264.

M S K **G** E E L F T G
 ATGTCTAA**GGT**GAAGAATTATTTACTGGT..




At the fourth position, the influence of amino acids was greater than that of synonymous codons.

Effect of N-terminal tag on the protein synthesis

PUREflex 2.1 (4 mM GSH)
 + template DNA
 incubation
 for 24 hours
 at 23, 30 or 37°C
 SDS-PAGE


Luc : MEDAKNIKKG..
 ALP: MRTPEMPVLE..
 DHFR: MISLIAALAV..

SKG/SKY



Fuse-Murakami et al. (2024) *Int. J. Mol. Sci.*, 25, 5264.

ALP




N-terminal additional sequence increased the amount of synthesized ALP at all temperature.

Expression of sfGFP variants in *E. coli*

BL21(DE3) + pET-sfGFP DNA
(Tyr/Asn/Ala/Gly@4th)
Incubation at 23, 30 or 37°C
+ IPTG (at OD=1)
Incubation at 23, 30 or 37°C
Measure OD and fluorescence

Fuse-Murakami *et al.* (2024) *Int. J. Mol. Sci.*, 25, 5264.



Variants that were highly synthesized in the PURE system were also highly expressed in *E. coli* cells.
The difference in the expression level between variants was less than that in the PURE system.
Expression of Ala and Gly variants suppressed the growth of *E. coli* cells.

Summary

Optimum sequence for the translation in the PURE system

- ▶ **AT-rich codons are preferable in the N-terminal region.**
- ▶ **Multiple codons should be used for each amino acid in the entire ORF.**
- ▶ **Amino acids at the N-terminal region are very important for efficient protein synthesis especially at lower temperature.**

Summary

Translation in the PURE system = Translation in *E. coli* cells ?

- ▶ Initiation factors are not essential.
- ▶ AT-rich region in the 5'UTR is important as well as SD sequence.
- ▶ Translation efficiency at lower temperature is dependent on the sequence of the target proteins.

Summary

Translation in the PURE system \approx Translation in *E. coli* cells

- ▶ Initiation factors are not essential.
- ▶ AT-rich region in the 5'UTR is important as well as SD sequence.
- ▶ Translation efficiency at lower temperature is dependent on the sequence of the target proteins.