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Exploration of optimal conditions for the incorporation of noncanonical amino acids by amber suppression using the PURE system

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
○Rena Matsumoto, Takayoshi Watanabe, Takahiro Hohsaka and Takashi Kanamori (1)GeneFrontier Corp., (2)School of Material Science, JAIST

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Incorporation of noncanonical amino acids by amber suppression using PUREfrefx[®]


Purpose of this study

Exploration of the optimal conditions for incorporation of a noncanonical amino acid into the N-terminal region by amber suppression

1. Optimal position of amber codon for incorporation of a noncanonical amino acid

Model protein: Stag-DHFR-FLAG

MNKL **KETAAAKFERSQHMS** DHFR FLAG



ver. 1 (ΔRF1)
+ 0.5 μg/μL amber suppressor tRNA aminoacylated with BODIPYFL-amino-phenylalanine (BFLAF-tRNA_{Cua})
1 ng/μL Stag(am)-DHFR-FLAG reaction at 37 °C for 60 min

The amount of the product containing BODIPY was higher at positions T7, A10 and F12, than at other positions. However, the efficiency of the incorporation of BFLAF was low.

2. Optimization of the concentration of three factors in the reagent for amber suppression

2-1. suppressor tRNA


ver. 1 (ΔRF1)
+ 0-2 μg/μL BFLAF-tRNA_{Cua}
1 ng/μL Stag(T7am or F12am)-DHFR-FLAG reaction at 37 °C for 60 min

The product containing BFLAF exhibited high expression levels when 0.5 or 1 μg/μL of suppressor tRNA was applied. However, the productivity was decreased when 2 μg/μL suppressor tRNA was used.

2-2. RF1

ver. 1 (ΔRF1)
+ 0/0.25 μg/μL BFLAF-tRNA_{Cua}
+ 0-10 μM RF1
1 ng/μL Stag(T7am or F12am)-DHFR-FLAG reaction at 37 °C for 60 min


The suppression efficiency was comparable to that observed in the absence of RF1, even at concentrations below 100 nM. Additionally, full-length products were synthesized without using suppressor tRNA.

The presence of RF1 is recommended to prevent the synthesis of full-length products incorporating near-cognate tRNA and to achieve precise amber suppression.

2-3. Magnesium

ver. 1 (ΔRF1/ΔCTP/ΔUTP/ΔT7RNAP)
+ 0.25 μg/μL BFLAF-tRNA_{Cua}
+ 0/0.1 μM RF1
Mg: 8-17 mM
0.2 μg Stag(T7am or A10)-DHFR-FLAG mRNA reaction at 37 °C for 60 min

The product containing BFLAF was synthesized with high efficiency at 14 mM magnesium or higher.

sup. tRNA: 0.5-1 μg/μL
RF1: 0-100 nM
Magnesium: > 14 mM

3. Optimal codon immediately before amber codon

In *E. coli* S30-extract system, it has been shown that the codon immediately preceding amber codon affects the efficiency of amber suppression using BFLAF.



Construction of DNA in which the codon immediately preceding amber codon was substituted

Stag (X6T7am)-DHFR-FLAG


MNKL **KX**AAAKFERSQHMS DHFR FLAG

X=F (ttt)
F (ttc)
L (tta)
L (ttg)
.....
*tag (amber codon)

PUREfrefx ver. 1 (ΔRF1)
+ 0.25 μg/μL BFLAF-tRNA_{Cua}
+ 100 nM RF1
1 ng/μL Stag(X6T7am)-DHFR-FLAG reaction at 37 °C for 60 min

$$RF_x (\%) = \frac{FL_x}{FL_{GAA-full} + FL_{GAA-short}} \times 100$$

The codon immediately upstream of amber codon significantly affected the synthesis by amber suppression using PUREfrefx.

Good codon: Ala (GCT/ GCC/ GCA/ GCG)
Gly (GGT/ GGC/ GGA/ GGG)
Arg (CGA)
Asn (AAT)

4. Conclusion

- A noncanonical amino acid was incorporated by amber suppression even in the presence of RF1 (or it could be better with RF1).
- The position of amber codon and a codon immediately before amber codon significantly affected the suppression efficiency.

For more information, please contact us.
URL: www.genefrontier.com
E-mail: purefrefx@genefrontier.com